

Indigo™ 202 Digital Transmitter

		1	2	3	4	5
Vaisala Indigo™ 202 Digital Transmitter		INDIGO202	A	0	A	0
1	Display	With metric units With non-metric units	1 2			
2	Serial protocol	Modbus RTU over RS-485¹	A			
3	For future use	None		0		
4	For future use	None			A	
5	Customization	None				0

Selections in bold are included in the prices of the basic version.

¹) Default settings: 19200 bps; databits:8; parity: None; stopbits: 2

Examples of order codes with typical settings:

1	2	3	4	5
INDIGO202	1	A	0	A 0

Accessories

Remote probe cable for Indigo,1 m
Remote probe cable for Indigo,3 m
Remote probe cable for Indigo,5 m
Remote probe cable for Indigo,10 m

Item code

INDIGOCABLE1M
INDIGOCABLE3M
INDIGOCABLE5M
INDIGOCABLE10M